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Human-Spongiosa, gefriergetrocknet, DIZG

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Fax: +49 (0)30 6576 3055

1. NAME OF THE MEDICINAL PRODUCT

Human cancellous bone, freeze-dried, DIZG

2. COMPOSITION

Human cancellous bone tissue.

3. PHARMACEUTICAL FORMS

Cubes, chips, granules, cylinders.

4. CLINICAL PARTICULARS

4.1 Therapeutic indications:

For implantation as tissue replacement for lost or insufficient tissue and for filling and stabilising bony defects. It is used in the field of general, neuro and accident surgery as well as in orthopaedics. Further applications in other operative disciplines are also conceivable.

4.2 Posology and method of administration:

For grafting in accordance with the clinical indication. The contents and outer surface of the inner bag are sterile. Prior to being used, human cancellous bone, freeze-dried, DIZG, is to be rehydrated at room temperature in an appropriate sterile, physiological medium for at least 30 minutes (e.g. in isotonic infusion solution).

4.3 Contraindications:

Grafting in necrotic operation sites is contra-indicated. Due to an expected drop in the settling rate, accurate diagnosis of grafting in undersupplied or infected operation sites is imperative.

4.4 Special warnings and precautions for use:

For single use only.

4.5 Interactions with other medicinal products and other forms of interaction:

None known.

4.6 Pregnancy and lactation:

There are no sufficient data on the use of human cancellous, freeze-dried, DIZG during pregnancy or lactation.

4.7 Effects on ability to drive and operate machines:

Not applicable.

4.8 Undesirable effects:

None known.

The following incidence rates are used to assess side effects: very common (>10%); common (>1% – <10%); uncommon (>0.1% – <1%); rare (>0.01% – <0.1%); very rare (<0.01% or unknown).

Going on case history and physical examination results, each donor is carefully checked for any reasons for exclusion. Extensive serological tests for infection and tests for any viruses are also performed using nucleic acid amplification techniques (NAT). Screening currently includes tests for antibodies to hepatitis A, B, C and HIV 1/2 as well as CMV, Treponema pallidum and HTLV I/II and NAT tests for hepatitis B, C and HIV genomes. These tests exceed the requirements of EU Directive 2006/17. Despite these extensive tests and the use of a validated chemical cold sterilisation process, the transmission of infectious diseases through previously unknown or unidentifiable pathogens cannot be ruled out with complete certainty.

4.9 Overdose:

Not applicable.

5. PHARMACOLOGICAL PROPERTIES

Pharmacotherapeutic group:

Graft of human origin.

5.1 Pharmacodynamic and 5.2 Pharmacokinetic properties:

Owing to the primarily physical mechanism of action thereof, studies investigating the pharmacodynamics and pharmacokinetics of this therapeutic agent in trial subjects or patients are not available. Positive results in respect of the effectiveness and tolerability of the preparation in 1205 patients have been substantiated in the „Study report on the clinical effectiveness and tolerability of allogenic, avital tissue implants sterilised in peracetic acid/ethanol“. This report is based on the statistically analysed and evaluated experience gained over a period of four years in the main hospitals using the preparation. It details excellent to good results in 90% of all cases in respect of the use of human cancellous as a spacer and restorative agent. Operating surgeons rated the settling as very good to good in 82% of all cases. Positive results (4 patients) were achieved in the bridging-over of defects in perioperative fractures among patients with status post implantation of titanium hip endoprostheses.

5.3 Preclinical safety data:

Given the decades of use in humans, experimental pre-clinical animal studies have not been conducted.

6. PHARMACEUTICAL PARTICULARS

6.1 List of excipients:

Not applicable.

6.2 Incompatibilities:

None known.

6.3 Shelf life:

5 years in the unopened container.

6.4 Special precautions for storage:

Do not store above 25 °C.

6.5 Nature and contents of container:

The contents and the external surface of the internal packaging are sterile.
cancellous cube:

1 piece 1 cm³, 3 cm³, 9 cm³, 12 cm³, 15 cm³

3 pieces 1 cm³, 3 cm³

cancellous chips <1 cm: 5 cm³, 10 cm³, 15 cm³, 20 cm³, 30 cm³

Cortical/cancellous chips <1 cm: 5 cm³, 10 cm³, 15 cm³, 30 cm³, 45 cm³, 60 cm³

Cortical/cancellous chip, J shape: 1 piece

10 x 10 – 15 mm, 4 – 5 mm thick;

15 x 20 mm, 10 mm thick

Granules:

Particle size 0.212 – 0.85 mm: 0.8 cm³, 1.6 cm³

Particle size 0.25 – 1 mm: 0.5 cm³, 1.0 cm³, 1.5 cm³,

Particle size 1 – 2 mm: 0.5 cm³, 1.0 cm³, 1.5 cm³,

cancellous cylinder: 1 piece

10 x 30 mm (height), 15 x 15 mm (height), 20 x 15 mm (height)

Caput femoris:

endoprothetic, halved, without cartilage: 1 piece and 2 pieces

endoprothetic, quartered, without cartilage: 1 piece and 2 pieces

halved, without cartilage: 1 piece and 2 pieces

quartered, without cartilage: 1 piece

Wedge 1 piece, 22.5°, 45°

Os ilium, bicortical: 1 piece

10 – 19 x 40 mm, 10 – 19 x 60 mm, 20 x 40 mm, 20 x 60 mm

Os ilium, tricortical: 1 piece

10 – 19 x 40 mm, 10 – 19 x 60 mm, 20 x 40 mm,

40 x 20 mm, 20 x 60 mm, 60 x 20 mm

Not all pack sizes may be marketed.

6.6 Special precautions for disposal and other handling:

The content of one pack is intended for use in one patient only. Any remaining tissue is to be discarded (disposal in accordance with AVV disposal key AS 180102).

Upon being opened, the graft must be used within 2 hours.

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10. DATE OF REVISION OF THE TEXT

02.07.2012

11. GENERAL CLASSIFICATION FOR SUPPLY

Medicinal product subject to medical prescription.

Adverse outcomes attributable to the tissue must be promptly reported to your local representative.

Documentation: patient file labels are enclosed with each graft for documentation purposes.

®

DIZG
German Institute for Cell and Tissue Replacement

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