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Demineralised human bone matrix, freeze-dried, DIZG

**2. COMPOSITION**

Demineralised human bone matrix, freeze-dried, DIZG, consisting of human decalcified bone tissue.  
 DBM pasty contains: 30% demineralised bone matrix (DBM) in a mixture of sodium hyaluronate, common salt, sodium monophosphate and sodium diphosphate in water for injection.  
 DBMx-press contains: 36% demineralised bone matrix (DBM) in a mixture of sodium hyaluronate, common salt, sodium monophosphate and sodium diphosphate in water for injection.

**3. PHARMACEUTICAL FORMS**

Chips, splinters, granules and paste.

**4. CLINICAL PARTICULARS****4.1 Therapeutic indications:**

For implantation. For the grafting of bony defects in the field of mouth and jaw surgery. Further applications in other operative disciplines are also conceivable.

**4.2 Posology and method of administration:**

For grafting in accordance with the clinical indication. The bone defects are to be filled loosely, as the bone matrix expands. Prior to being used, demineralised human bone matrix, freeze-dried, DIZG, is to be rehydrated at room temperature in an appropriate sterile, physiological medium for at least 30 minutes (e.g. in isotonic infusion solution).

**Note:** DBM pasty ad DBMx-press do not require any rehydration!

**4.3 Contraindications:**

Grafting in necrotic operation sites is contra-indicated. Due to an expected drop in the settling rate, accurate diagnosis of grafting in undersupplied or infected operation sites is imperative.

**4.4 Special warnings and precautions for use:**

For single use only.

**4.5 Interactions with other medicinal products and other forms of interaction:**

None known.

**4.6 PREGNANCY AND LACTATION:**

There are no sufficient data on the use of demineralised human bone matrix, freeze-dried, DIZG during pregnancy or lactation.

**4.7 EFFECTS ON ABILITY TO DRIVE AND OPERATE MACHINES:**

Not applicable.

**4.8 UNDESIRABLE EFFECTS:**

None known.

The following incidence rates are used to assess side effects: very common (>10%); common (>1% – <10%); uncommon (>0.1% – <1%); rare (>0.01% – <0.1%); very rare (<0.01% or unknown).

Going on case history and physical examination results, each donor is carefully checked for any reasons for exclusion. Extensive serological tests for infection and tests for any viruses are also performed using nucleic acid amplification techniques (NAT). Screening currently includes tests for antibodies to hepatitis A, B, C and HIV 1/2 as well as CMV, Treponema pallidum and HTLV I/II and NAT tests for hepatitis B, C and HIV genomes. These tests exceed the requirements of EU Directive 2006/17. Despite these extensive tests and the use of a validated chemical cold sterilisation process, the transmission of infectious diseases through previously unknown or unidentified pathogens cannot be ruled out with complete certainty.

**4.9 Overdose:**

Not applicable.

**5. PHARMACOLOGICAL PROPERTIES****Pharmacotherapeutic group:**

Graft of human origin.

**5.1 Pharmacodynamic and 5.2 Pharmacokinetic properties:**

Owing to the primarily physical mechanism of action thereof, studies investigating the pharmacodynamics and pharmacokinetics of this therapeutic agent in trial subjects or patients are not available. Positive results in respect of the effectiveness and tolerability of the preparation in 36 patients have been substantiated in the „Study report on the clinical effectiveness and tolerability of allogenic, avital tissue implants sterilised in peracetic acid/ethanol“. This report is based on the statistically analysed and evaluated experience gained over a period of four years in the main hospitals using the preparation. It details excellent to good results throughout in respect of the use of human DBM together with autologous bone tissue and platelet-rich plasma for the augmentation of upper and lower jaw defects. Reference is also

made to a placebo-controlled study, in which excellent defect grafting and excellent settling were able to be achieved in 50 patients.

**5.3 Preclinical safety data:**

Given the decades of use in humans, experimental pre-clinical animal studies have not been conducted.

**6. PHARMACEUTICAL PARTICULARS****6.1 List of excipients:**

DBM pasty contains: 30 % demineralised bone matrix (DBM) in a mixture of sodium hyaluronate, common salt, sodium monophosphate and sodium diphosphate in water for injection.

DBMx-press contains: 36 % demineralised bone matrix (DBM) in a mixture of sodium hyaluronate, common salt, sodium monophosphate and sodium diphosphate in water for injection.

**6.2 INCOMPATIBILITIES:**

None known.

**6.3 SHELF LIFE:**

5 years in the unopened container.

**6.4 SPECIAL PRECAUTIONS FOR STORAGE:**

Do not store above 25 °C.

**6.5 NATURE AND CONTENTS OF CONTAINER:**

The contents and the external surface of the internal packaging are sterile.

**Granules:**

**Particle size 0.08 – 0.8 mm:** 1 cm<sup>3</sup>, 1.5 cm<sup>3</sup>, 3 cm<sup>3</sup>, 5 cm<sup>3</sup>, 10 cm<sup>3</sup>, 30 cm<sup>3</sup>

**Particle size 0.212 – 0.850 mm:** 1 cm<sup>3</sup>, 2 cm<sup>3</sup>, 3 cm<sup>3</sup>

**DBM in DBM putty:**

0.15 cm<sup>3</sup> in 0.5 cm<sup>3</sup> DBM putty

0.3 cm<sup>3</sup> in 1.0 cm<sup>3</sup> DBM putty

0.75 cm<sup>3</sup> in 2.5 cm<sup>3</sup> DBM putty

1.5 cm<sup>3</sup> in 5.0 cm<sup>3</sup> DBM putty

3 cm<sup>3</sup> in 10.0 cm<sup>3</sup> DBM putty

**DBM in DBMx-press:**

0.9 cm<sup>3</sup> (1 – 6 mm) in 2.5 cm<sup>3</sup> DBMx-press

1.8 cm<sup>3</sup> (1 – 6 mm) in 5 cm<sup>3</sup> DBMx-press

3.6 cm<sup>3</sup> (1 – 6 mm) in 10 cm<sup>3</sup> DBMx-press

5.4 cm<sup>3</sup> (1 – 6 mm) in 15 cm<sup>3</sup> DBMx-press

7.2 cm<sup>3</sup> (1 – 6 mm) in 20 cm<sup>3</sup> DBMx-press

**Particle size 1 – 2 mm:** 1 cm<sup>3</sup>, 1.5 cm<sup>3</sup>, 3 cm<sup>3</sup>, 5 cm<sup>3</sup>

**Particle size 1 – 3 mm:** 1 cm<sup>3</sup>, 1.5 cm<sup>3</sup>, 3 cm<sup>3</sup>, 5 cm<sup>3</sup>, 10 cm<sup>3</sup>, 15 cm<sup>3</sup>, 30 cm<sup>3</sup>

**Particle size 4 – 5 mm:** 5 cm<sup>3</sup>, 10 cm<sup>3</sup>, 15 cm<sup>3</sup>, 20 cm<sup>3</sup>, 30 cm<sup>3</sup>

**Particle size 6 – 10 mm:** 5 cm<sup>3</sup>, 10 cm<sup>3</sup>, 15 cm<sup>3</sup>, 30 cm<sup>3</sup>

**Demineralised cortical strut, flex, 2 – 3 mm thick:**

15 x 25 mm, 25 x 25 mm, 15 x 50 mm, 15 x 100 mm, 15 x 150 mm

Not all pack sizes may be marketed.

**6.6 Special precautions for disposal and other handling:**

The content of one pack is intended for use in one patient only. Any remaining tissue is to be discarded (disposal in accordance with AVV disposal key AS 180102). Upon being opened, the graft must be used within 2 hours.

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**10. DATE OF REVISION OF THE TEXT**

07.02.2013

**11. GENERAL CLASSIFICATION FOR SUPPLY**

Medicinal product subject to medical prescription.

Adverse outcomes attributable to the tissue must be promptly reported to your local representative

**Documentation:** patient file labels are enclosed with each graft for documentation purposes.

**DIZG**  
**German Institute for Cell and Tissue Replacement**  
 Koepenicker Strasse 325  
 D-12555 Berlin, Germany  
 Phone: +49 (0 30) 65 76 30 50  
 Fax: +49 (0 30) 65 76 30 55

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